## COMPOSITION

Hydroxyde de calcium, charge hydrophile, plastifiant

## PROPRIETES

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EN

**COMPOSITION**

Calcium hydroxide, hydrophilic filler, plasticizer

**PROPERTIES**

Calcium hydroxide forms an alkaline medium (pH more than 12) which provides long antibacterial effect and results in formation of dental-cement barrier. That creates conditions for growth and formation of tooth root, regeneration of bone and periapical tissues, increases pulp regeneration ability and stimulates formation of secondary dentine.

Endodontic paste **CALCETIN ENDO**, for temporary filling of root canals, is produced as a radiopaque filler, which reduces perceptible changes within 3–5 weeks, if changing the paste 2–4 times during treatment.

There are no complications in case of accidental overfilling of tooth canal. Radiopaque filler allows to control the material in the canal at X-ray.

**INDICATIONS**

Paste **CALCETIN ENDO** is used for:

- treatment of pulpitis in permanent and temporary teeth with incomplete roots growth and formation (apexogenesis)
- formation of osteoid-cement barrier in not shaped top or root of the tooth during:
  - treatment of irreversible forms of pulpitis
  - conservative treatment of all forms of chronic periodontitis
  - in case of any root perforation in permanent teeth with shaped roots
- treatment of canals with not shaped apex in temporary teeth at periodontitis

**CONTRAINDICATIONS**

Hypersensitivity to components of medical device.

Do not misuse.

**SIDE EFFECTS**

As far as all terms of storing, transporting and application are observed, there are no side effects.

**METHOD OF USE**

**Attention!** It is recommended to store material at room temperature (22–24) °C within 2 hours before use.

During treatment of pulpitis with not shaped roots in permanent teeth, apply paste **CALCETIN ENDO** directly on the tooth pulp. Isolate tooth cavity with temporary strong cement. In 3–4 weeks change the paste for fresh portion. For further treatment, paste should be changed in accordance with X-ray and EOD results every 2–3 month till complete growth and formation of the root. After tooth root top formation, fill root canal with permanent cement and restore the tooth crown.

For apexification of teeth with uncompleted root growth, tool and drug treatment of root canal should be fulfilled. Insert paste **CALCETIN ENDO** for endodontic treatment in dried canal and isolate with temporary strong cement. Treatment is carried out under X-ray control. In 3–4 weeks it is recommended to change the paste in the canal for fresh portion. After 2–3 month of treatment check parodontal condition with parodontal explorer or X-ray. Apply fresh portion of paste once per 3 month till complete formation of apical barrier. After complete apical barrier formation, fill the canal with sealer.

For treatment of chronic periodontitis, fulfill tool and drug treatment of the canal, after diagnostic testing of morphological features of tooth roots and periapical tissues (with the help of X-ray). Insert paste **CALCETIN ENDO** for endodontic use in prepared tooth canal. Isolate the paste with dry cotton ball and apply temporary filling for 7–10 days. Depending on the results, continue treatment procedures with fresh portion of the paste or fill the canal with permanent sealer.

During conservative treatment of chronic periodontitis and any perforations of root canals after application of paste for endodontic treatment **CALCETIN ENDO** for 3–6 weeks, it is possible to fill the canal with permanent sealer.

**Attention!** Wipe the tip of a syringe with cotton ball and close syringe cap after use, since paste hardens in syringe tip at contact with air.

**PRESENTATION**

Paste (syringe)	2 g / 2 x 2 g
Tips	2 pc or 6 pc

**STORAGE**

Store at the temperature from 5 °C to 25 °C.

Do not use after expiry date.

Shelf life — 3 years

**PRECAUTIONS**

Wear suitable protective clothing, gloves and goggles before use.

Avoid contact with eyes. In case of eye or skin contact, as well as contact with soft tissues and mucosa of oral cavity, remove the material immediately and wash with water. Use with caution in patients with mental disorders.

Do not use material for pregnant women (because of insufficient clinical data).

In cases of allergic reactions in particularly sensitive patients, the material should be removed and its further application abandoned.

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Tehnodent Co. Ltd,