

# INSTRUCTION FOR USE CRESOTIN

## Paste

### COMPOSITION

Zinc oxide, zinc sulfate, radiopaque filler — barium sulfate, dexamethasone, paste-forming agent

### PROPERTIES

Zinc oxide and zinc sulfate are antiseptics, having astringent properties. Barium sulfate is a radiopaque filler. Antiseptic base of the paste — chlorophenol and camphor — allows applying paste after treatment with CRESOTIN liquids (TehnoDent) as a finishing part of treatment with materials of the same antiseptic base. Paste contains dexamethasone (0.01 %) for reducing pain when inserting material out of the apex.

Paste hardens in dentinal medium (in-the-canal moisture).

Paste has an excellent fluidity (23 mm), easily inserted into canal with the help of root canal filling instrument and hardens under the influence of in-the-canal moisture during 72 hours. After hardening, the solubility of a material is 1.5 %.

### INDICATIONS

Paste **CRESOTIN** is used for filling of infected and complex root canals as well as for filling canals with incomplete pulp extirpation.

### CONTRAINDICATIONS

Drug idiosyncrasy. Do not misuse.

### SIDE EFFECTS

As far as all terms of storing, transporting and application are observed, there are no side effects.

### METHOD OF USE

**Attention!** It is recommended to store material at room temperature (22–24) °C within 1–2 hours before use. During treatment of multi-canal teeth, dry thoroughly the pulp cavity and canals before filling.

After thorough mechanical treatment of root canal, fill it with **CRESOTIN** paste with the help of root canal filling instrument. When filling canal with the use of gutta-percha points, first apply the paste on gutta-percha point and after insert it in the canal.

It is not recommended to insert material out of root apex.

For quicker hardening of the paste, canal can be treated with distilled water-wet paper point.

After filling of root canal, it is necessary to isolate it with GI cement and make a further restoration of the crown.

**Warning!** In case of allergic reactions with hypersensitive patients, it is necessary to extract the paste from the canal and stop further use of material.

**PRESENTATION**

Paste (jar)

25 g

**STORAGE**

Store at the temperature from 5 °C to 25 °C. Keep in dry place.

Tightly close container immediately after use.

Do not use after expiry date.

Shelf life — 3 years.

**MANUFACTURER**

TehnoDent Co. Ltd

1/6 Beryozovaya str., Severnyi,  
Belgorod region, Russia, 308570

Tel. +7 (495) 118-28-73

E-mail: [info@tehnodent.org](mailto:info@tehnodent.org)

[www.tehnodent.org](http://www.tehnodent.org)

***For professional use only***

1.015.02-03\_EN\_002

Rev.: 2019-09