



RELIGHT *colour*

LIGHT-CURING COLOUR

The product is designed to imitate colour effects (chalky white spots, enamel cracks, etc.) and tone the composite restoration

PASTE (1 g)

*Designed
to Inspire...*

TehnoDent company, 2022

+7 495 118-28-73
www.tehnodent.org
info@tehnodent.org

- White
 - Chalk spots imitation
 - Dental cusps and rims emphasis
 - Covering of pigmented tissues
- Yellow
 - Emphasizing of fissures and cervical area
 - Toning of dentine
- Brown
 - Imitation of cracks and pigmented fissures
- Grey
 - Imitation of transparent enamel area
- Ochre
 - Emphasis of fissures and cervical area